

Endogenous GPCR Transcript Profiling of Life Technologies' Primary Cell Offerings: Opening the Door to Physiologically Relevant Cell Systems for Pharmaceutical Screens



Rhonda A. Newman & David T. Kuninger, Life Technologies, Primary & Stem Cells Business Unit, 29851 Willow Creek Rd, Eugene, OR, USA 97402

ABSTRACT

The association of ligands with cell surface G protein-coupled receptors (GPCRs) elicits a myriad of physiological responses through activation of intracellular signaling cascades. Mutations in GPCRs and associated G-proteins have been implicated in a number of diseased states. Given the localization of these receptors and their role in pathogenesis, it is not surprising that approximately fifty percent of all pharmaceuticals target this family of receptors.

Heterologous over-expression of GPCRs is commonly utilized to screen the efficacy of GPCR-targeted drug candidates. While these engineered cell lines have certain utility, several significant problems have been noted: receptors (1) are expressed at levels far-exceeding physiologically relevant concentrations, (2) can hetero-oligomerize and/or cross talk with non relevant endogenous receptors, and (3) many times are not appropriately targeted to the plasma membrane. In addition, intellectual property places restrictions on the use of certain GPCR expressing cell lines. Together, these caveats underscore the need for alternative cellular assay systems that overcome these limitations in the discovery and development of GPCR-targeted pharmaceuticals.

Primary cells provide physiologically relevant systems for the assessment of complex signaling events, with appropriately targeted endogenous receptors and without the restrictions of IP. Yet, use of primary cells in GPCR drug development has been limited due to several factors that can include relatively low endogenous receptor expression and/or presence of multiple related receptor isoforms that may complicate analysis. Additionally, donor-to-donor variability and issues of scale are important considerations. To address these issues, we undertook GPCR transcript profiling of normal neonatal human epidermal keratinocytes (HEKn) using TaqMan® arrays and quantified GPCR transcripts levels across a range of donors and through multiple passages.

The GPCR TaqMan® array available from Life Technologies/Applied Biosystems targets 367 GPCR transcripts representing 50 subfamilies. Our results show remarkable consistency in transcript levels both across HEKns from 2 donors and during expansion throughout 6 passages (18 population doublings). Among the subfamilies present within our primary keratinocytes are LPA receptors, adrenoreceptors, retinoic acid receptors, and proteinase-activated receptors. Many of these endogenous receptors have been, and continue to be, targets of pharmaceutical development to treat tumor metastasis, pain, and neurological and inflammatory disorders. Of particular interest are orphan receptors that are expressed within our primary keratinocytes, including receptors which have been implicated in tumorigenesis (GPR56 and GPR110). While the endogenous ligands of these receptors have yet to be identified, they represent interesting targets for basic research and pharma development. Taken together our data highlight the value of using physiologically relevant primary human cell models in drug discovery and development.

MATERIALS AND METHODS

Total RNA was isolated from primary cells grown in varying basal medium formulations available from Life Technologies. Using cDNA libraries, created by reverse transcription, the GPCR transcript levels were determined via TaqMan® array analysis and it was determined whether or not GPCR transcript levels change as a function of donor-to-donor variation, upon increasing passage, between epithelial cell types grown in various basal media and supplements, and between epithelial and mesenchymal cell types.


RESULTS

Figure 1. GPCR Transcript Levels Within 4' Neonatal Keratinocytes Were Assessed using the TaqMan® Gene Signature Array for Human GPCRs

GPCR	HEKn (Epilife + EDGS)	HCEC (Epilife + HCGS)	HCEC (KSFM)	HDFn (DMEM + 10% FBS)
ADORA2B	30.7	31.7	32.4	
ADRB2	29.2	31.3	32.2	33.8
BAI2				34.7
BDRB1	32.5	34.6		32.4
GPR137	31.4	32.3	33.1	31.8
CALCR	34.8			28.9
CCRL1				29.9
CB1R	32.6	32.6	32.5	29.9
CELSR1	30.6	32.8	32.7	
CELSR2	33.0			33.4
CHRM2				33.9
DRD5	34.4			33.9
EPOR				32.5
LPAR1	30.6	31.5	32.2	29.9
LANCL1				30.4
OPN3				31.7
SIPR2	32.9			32.6
TRPV1	29.9	32.5	32.6	
SIPR5	31.3	31.6		
EDNRB	34.8		34.3	34.8
EMR2	34.9			33.9
ELTD1				33.9
FZD1	29.5	32.0	32.6	31.4
FZD10	34.0			33.0
FZD16	34.2			33.0
FZD2	32.9	35.0		31.8
FZD3	33.2	33.2	33.2	
FZD4				31.6
FZD5	31.8	32.8	33.9	34.8
FZD6	29.4	30.4	30.6	31.2
FZD7	33.0	34.7		29.7
GABBR1	31.1	33.9	33.4	32.9
GPR106	29.5	31.5	32.6	28.7
GPR1				31.1
GPR110	31.7	32.8	33.1	
GPR115	32.6	34.1	34.9	
GPR124				29.6
GPR125	31.8	33.3	33.8	34.8
GPR126	31.3	33.9	34.2	
GPR135	34.8			
GPR153	29.2	30.6	30.8	31.4
GPR161	32.0	33.9	33.9	32.5
GPR181				32.8
TRPV5	32.3			



- GPCR transcript levels were assessed for 367 GPCRs from 50 subfamilies
- 30 GPCRs (~10% of profiled GPCRs) were shown to have cycle thresholds (C_T) ≤ 32
- Among the subfamilies represented were LPA receptors, adrenoreceptors, retinoic acid receptors, and proteinase-activated receptors
- Many orphan receptors were also identified including receptors that have been implicated in tumorigenesis (GPR56 and GPR110)

Figure 2. Keratinocytes Can Be Passaged Multiple Times with Few Significant Changes in GPCR Transcript Levels



- Only 2 GPCRs were shown to have > 2-fold changes in transcript levels, allowing large expansion of primary banks of cells for use in primary HTS screens.

Figure 3. GPCR Transcript Levels Within Keratinocytes Vary Little from Donor to Donor



- No additional GPCRs were identified in the subsequent screening of an additional lot of keratinocytes using the Taqman® Gene Signature Array for Human GPCRs.
- Only 5 GPCRs were shown to have > 2-fold changes in transcript levels.

Figure 4. Other Primary Cell Offerings Available From Life Technologies Preliminarily Screened for GPCR Transcripts




Figure 5. Preliminary GPCR Transcript Analysis of Epithelial (HEKn & HCEC) and Mesenchymal (HDFn) Cell Types Available From Life Technologies

GPCR	HEKn (Epilife + EDGS)	HCEC (Epilife + HCGS)	HCEC (KSFM)	HDFn (DMEM + 10% FBS)
ADORA2B	30.7	31.7	32.4	
ADRB2	29.2	31.3	32.2	33.8
BAI2				34.7
BDRB1	32.5	34.6		32.4
GPR137	31.4	32.3	33.1	31.8
CALCR	34.8			28.9
CCRL1				29.9
CB1R	32.6	32.6	32.5	29.9
CELSR1	30.6	32.8	32.7	
CELSR2	33.0			33.4
CHRM2				33.9
DRD5	34.4			33.9
EPOR				32.5
LPAR1	30.6	31.5	32.2	29.9
LANCL1				30.4
OPN3				31.7
SIPR2	32.9			32.6
TRPV1	29.9	32.5	32.6	
SIPR5	31.3	31.6		
EDNRB	34.8		34.3	34.8
EMR2	34.9			33.9
ELTD1				33.9
FZD1	29.5	32.0	32.6	31.4
FZD10	34.0			33.0
FZD16	34.2			33.0
FZD2	32.9	35.0		31.8
FZD3	33.2	33.2	33.2	
FZD4				31.6
FZD5	31.8	32.8	33.9	34.8
FZD6	29.4	30.4	30.6	31.2
FZD7	33.0	34.7		29.7
GABBR1	31.1	33.9	33.4	32.9
GPR106	29.5	31.5	32.6	28.7
GPR1				31.1
GPR110	31.7	32.8	33.1	
GPR115	32.6	34.1	34.9	
GPR124				29.6
GPR125	31.8	33.3	33.8	34.8
GPR126	31.3	33.9	34.2	
GPR135	34.8			
GPR153	29.2	30.6	30.8	31.4
GPR161	32.0	33.9	33.9	32.5
GPR181				32.8
TRPV5	32.3			

- Analysis of transcript levels from epithelial cell types indicate significant changes (i.e., >> 2-fold) in GPCR transcripts between HEKns and HCECs.
- Analysis of transcript levels from HCECs grown in varying medium conditions indicate little to no change in GPCR transcript levels.
- Analysis of transcript levels within epithelial vs. mesenchymal cell types indicates significant differences in endogenous transcript levels.

CONCLUSIONS

- Many of the endogenous receptors identified in this study have been, and continue to be, targets of pharmaceutical development to treat tumor metastasis, pain, neurological and inflammatory disorders.
- The transcript levels of HEKn endogenous receptors are shown to be surprisingly consistent upon increasing passage of the cells, as well as from donor-to-donor, demonstrating the ability to expand primary cell banks to accommodate the scale required for primary HTS screens.
- Preliminary assessment of the GPCR transcripts of other primary cell offerings available from Life Technologies demonstrates significant differences in GPCR expression between epithelial and mesenchymal cell types.

ACKNOWLEDGEMENTS

We would like to thank Keith Olson, Trillium Blackmer, Bonnie Hanson, Anne Dickson, Mike O'Grady, and Trey Knollman for helpful discussions.