

NA-Fluor™ and NA-XTD™ Influenza Neuraminidase Assays for Neuraminidase Quantitation and Inhibition Assays



Anthony C. Chiulli, Corinne E. Miller, Albana Mihali, Laura Thibodeau, Rouh-Rong Juo, Brian D'Eon, and Melissa A. Gee
Life Technologies Corporation, Molecular & Cell Biology Division, 35 Wiggins Avenue, Bedford, MA 01730, USA

OVERVIEW

Applied Biosystems® NA-Fluor™ and NA-XTD™ Influenza Neuraminidase Assays:

- Applications:
 - Measure virus sensitivity to neuraminidase inhibitors (NIs)
 - Identify new candidate NIs or anti-viral compounds
- Results:
 - High sensitivity and wide dynamic range
 - Neuraminidase inhibitor IC₅₀ values in benchmarking to traditional assay protocols
 - Signal kinetics and stability provides assay format flexibility and ease-of-use
 - Excellent Z' values for high throughput screening assay capability
 - Cell-based viral growth inhibition assay demonstrated

INTRODUCTION

The Applied Biosystems® NA-Fluor™ and NA-XTD™ Influenza Neuraminidase Assay Kits are optimized for direct quantitation of neuraminidase inhibitor (NI) resistance by monitoring neuraminidase (NA) activity. These assays can also be used in the identification or characterization of novel NIs or other antiviral compounds during lead discovery. Global monitoring of neuraminidase activity is essential for understanding NI efficacy for seasonal, pandemic or avian influenza, and studying the epidemiology of resistance mutations. The rapid spread of seasonal H1N1 virus oseltamivir-resistance to 98.5% in the US (1) highlights the importance of resistance monitoring and the need for continued development of new influenza anti-virals.

The fluorescent NA-Fluor™ assay uses the traditional 2-[4-methylumbelliferyl]-D-N-acetylneuraminic acid (MUNANA) substrate. Assay reagents and protocol were optimized based on Neuraminidase Inhibitor Susceptibility Network (NISN) protocols to provide an economical, standardized and easy-to-use assay kit that generates IC₅₀ data comparable to historical MUNANA assay data. The NA-Fluor™ assay signal is stable up to several days after assay termination, enabling real-time flexibility and provides a typical Z' value of 0.8, for excellent high throughput screening assay capability. In addition, the NA-Fluor™ assay can be conducted in a real-time format for enzyme kinetic studies.

The chemiluminescent NA-XTD™ assay provides the next-generation NA-XTD™ 1,2-dioxetane neuraminidase substrate. Like the first-generation chemiluminescent NA-Star® Influenza Neuraminidase Inhibitor Resistance Detection Kit, the NA-XTD™ assay provides highly sensitive neuraminidase quantitation with a wide assay dynamic range. The NA-XTD™ assay provides extended-glow light emission. IC₅₀ values determined either immediately or several hours after assay completion are identical, and Z' values of 0.5 – 0.86 are obtained, providing excellent high throughput screening assay capability. The NA-XTD™ assay is also used to quantitate influenza NA activity directly in cell-based virus cultures to monitor viral growth or inhibition.

In addition to NI resistance quantitation (1-7), additional applications include identification of new NI compounds (8), NI characterization (9), studies of virus transmission (10), neuraminidase-based quantitation of virus-like particles (11), and cell-based virus quantitation (12).

MATERIALS AND METHODS

The NA-Fluor™, NA-XTD™ and NA-Star® assay kits are supplied by Life Technologies. The NA-Fluor™ assay (PN 4457091) includes 1) NA-Fluor™ 2X Assay Buffer, 2) NA-Fluor™ Substrate (solid powder) and 3) NA-Fluor™ Stop Solution. The NA-XTD™ assay (PN 4457635) includes 1) NA Sample Prep Buffer (contains Triton X-100), 2) NA-XTD™ Assay Buffer, 3) NA-XTD™ Substrate concentrate, 4) NA-XTD™ Accelerator and 5) NA-Star™ Detection Microplates (solid-white 96-well microplate). Influenza strains include type A/H1N1 (VR-1469™ (human), VR-1520™ (human), VR-1882™ (swine)) and B/Leu40 (VR-1535) from ATCC (Manassas, VA) and A/H1N1/Texas/36/91, wild-type sensitive and H275Y mutant oseltamivir-resistant strains, (kindly provided by the CDC, Atlanta GA). Assays were performed according to their respective protocols.

NA-Fluor™ Assay: 25 µL diluted virus (in 1X Assay Buffer)
25 µL NI dilution (in 1X Assay Buffer), pre-incubate x 30 min at 37°C
50 µL NA-Fluor™ substrate, incubate x 60 min at 37°C
100 µL Stop Solution, read with Ex 355 nm/Em 460 nm

NA-XTD™ Assay: 25 µL diluted virus (in Assay Buffer)
25 µL NI dilution (in Assay Buffer), pre-incubate x 30 min at 37°C
25 µL 1:1000 NA-XTD™ substrate (in Assay Buffer), incubate x 30 min at r.t.
60 µL Accelerator, read 1 sec/well

Fluorescent assays and chemiluminescent were both performed with the SpectraMax M5 multiplate microplate reader (Molecular Devices, Sunnyvale, CA). IC₅₀ analysis was performed using GraphPad Prism™ dose-response curve-fitting.

ACKNOWLEDGEMENTS

Influenza A/Texas/36/91 (H1N1) wild-type and H275Y strains were kindly provided by Dr. Larisa Gubareva (Influenza Branch, CDC, USA).
Oseltamivir carboxylate was kindly provided by F. Hoffmann-La Roche Ltd (Basel Switzerland), and zanamivir was kindly provided by GlaxoSmithKline (Research Triangle Park, NC).

RESULTS

Figure 1. Neuraminidase Inhibitor IC₅₀ Determination Assay Workflow






Figure 2. Sensitivity Comparison of NA-Fluor™, NA-XTD™ and NA-Star® Assays




Assay workflow for use of either the NA-Fluor™ or NA-XTD™ assay for NI IC₅₀ determination or lead identification with viral isolates.

Figure 3. Comparison of NA-Fluor™ Assay to Other Assays for IC₅₀ Determination



Protocols were run in parallel to determine IC₅₀ values of influenza A/Texas/36/91 (H1N1) oseltamivir-sensitive and resistant H275Y strains for oseltamivir carboxylate. The NA-Fluor™ assay provides IC₅₀ values and sensitivity similar to NISN-published, MUNANA-based protocols (12) and can be compared with historical data sets based on these protocols. IC₅₀ values obtained with the NA-XTD™ and NA-Star® assays are similar, but can be slightly lower than values obtained with MUNANA-based assays.

Figure 4. NA-Fluor™ Assay: Signal Stability Provides Read-Time Flexibility




NA-Fluor™ signal read-out remains stable for hours and up to days providing flexibility and batch-mode capability for high throughput applications. A) Signal stability at room temperature (25°C) following addition of NA-Fluor™ Stop Solution. B) Signal stability with plate storage at 4°C following stop solution addition. Fluorescent signal shown as % of the immediate assay read-out. IC₅₀ values determined using data collected up to 3 days after assay completion (plates stored at 4°C) are identical (data not shown).

Table 1. NA-Fluor™ Assay: Suitability for High Throughput Compound Screening

Virus Strain/Dilution	Read Time	Signal/Noise (NI RFU/ANI RFU)	%CV (±NI)	%CV (NI)	Z'	HTS Suitability
VR-1520/1:10	Immediate	8.5	4.3	5.4	0.78	Excellent
	2 hr	8.4	3.7	5.1	0.7	Excellent
	3 hr	8.4	3.6	4.9	0.8	Excellent
	4 hr	8.4	3.6	4.9	0.8	Excellent

The NA-Fluor™ assay demonstrates low intra-assay signal variability and provides a Z' value of 0.78 or greater, consistent across multiple read times, indicating that it has excellent quality as a high throughput screening assay for identification of candidate neuraminidase inhibitor compounds.

Figure 5. NA-Fluor™ Assay: Real-Time vs. End-Point Assay Format



The NA-Fluor™ assay can be run in either end-point or real-time kinetic mode. A) Substrate turn-over rate is linear for >2 hours enabling flexibility in assay format. Real-time acquired RFU values are typically 5-6 fold lower than RFUs acquired after addition of stop solution at the same time point. B) Oseltamivir IC₅₀ values determined with either real-time data by slope analysis (no stop solution) or a 60 minute end-point read-out (with stop solution) read-time for influenza AWS/33 (H1N1) or influenza B/Leu40 are nearly identical.

Figure 6. NA-Fluor™ Assay: Compatibility with Viral Inactivation Methods






Figure 7. NA-XTD™ Assay: Extended-Glow Light Emission Kinetics



The half-life of light emission following addition of NA-XTD™ Accelerator (TA) with the NA-XTD™ assay is ~2 hours, compared to ~10 minutes with the NA-Fluor™ assay. The longer-lived light emission kinetics eliminates the need for accelerator injection and luminometer instrumentation equipped with on-board injectors.

Figure 8. NA-XTD™ Assay: Read-Time Flexibility




IC₅₀ values were determined using data collected up to 5 hours after addition of NA-XTD™ Accelerator. The IC₅₀ curves and values are nearly identical at each time point, using influenza B/Leu40.

Table 2. NA-XTD™ Assay: Suitability for High Throughput Compound Screening

Virus Strain/Dilution	Read Time (min)	Signal/Noise (NI RFU/ANI RFU)	%CV (±NI)	%CV (NI)	Z'	HTS Suitability
VR-1520/1:10	100	36	13.9	4.6	0.86	Excellent
	200	36	13.9	4.6	0.86	Excellent
	300	36	13.9	4.6	0.86	Excellent
	400	36	13.9	4.6	0.86	Excellent
	500	36	13.9	4.6	0.86	Excellent
	600	36	13.9	4.6	0.86	Excellent
	700	36	13.9	4.6	0.86	Excellent
	800	36	13.9	4.6	0.86	Excellent
	900	36	13.9	4.6	0.86	Excellent
	1000	36	13.9	4.6	0.86	Excellent


The NA-XTD™ assay provides Z' values of 0.5-0.86, determined with two different virus strains at multiple virus concentrations, and similar across multiple read times, indicating excellent quality as a high throughput screening assay for identification of candidate neuraminidase inhibitor compounds.

Figure 9. NA-XTD™ Assay: Compatibility with Triton® X-100 Addition to Virus Sample



Triton X-100 detergent at 1% has been shown to inactivate flu virus while increasing neuraminidase activity (13). A) NA Sample Prep Buffer (containing 10% Triton X-100) was added to virus stocks (0.1% final), virus was diluted, and IC₅₀ values determined with oseltamivir carboxylate. IC₅₀ values are very similar following addition of Triton X-100. B) NA activity of virus dilutions + Triton X-100. Typically, a 2-4-fold increase in NA activity is obtained with Triton X-100 addition to virus stock.

Figure 10. NA-XTD™ Assay: Cell-Based Virus Growth Inhibition Assay



The NA-XTD™ Assay was used to monitor cell-based virus growth inhibition similar to the described AVINA (Accelerated Viral Inhibition with NA as read-out) assay (12). MDCK 96-well microplate cell cultures were infected with indicated virus strains, followed by incubation in the presence of varying concentrations of oseltamivir carboxylate. Samples of culture media were assayed 24 hours later. Quantitation of NA activity with the NA-XTD™ assay demonstrates inhibition of viral growth by oseltamivir in cell culture. This cell-based assay capability can be used for identification of compounds inhibiting virus growth/replication or quantitation of antibody-based viral growth inhibition.

CONCLUSIONS

The Applied Biosystems® NA-Fluor™ Influenza Neuraminidase Assay Kit, based on the widely-used fluorescent MUNANA substrate, provides a standardized, easy-to-use, cost-effective assay system ideal for high throughput or compound screening. The NA-XTD™ Influenza Neuraminidase Assay Kit is a next-generation chemiluminescent neuraminidase assay providing high assay sensitivity and 'glow' light emission kinetics for improved ease-of-use and high throughput capability.

Together these assays offer:

- Standardized reagents and protocols for:
 - Neuraminidase inhibitor (NI) resistance screening
 - High throughput identification of new NIs
- Choice of fluorescent or chemiluminescent detection technology
- Simple instrumentation requirements
- High sensitivity for use with low virus concentrations
- Compatibility with batch-mode processing and large-scale assay throughput required for HTS
- Excellent assay quality (Z') for high-throughput screening new NI lead identification
- Compatibility with virus inactivation reagents
- Broad specificity of influenza detection
- Flexibility in assay format and read-time
- Additional NA assay applications
 - Cell-based viral growth and inhibition assays
 - Detection of NA from other organisms (bacteria, viruses)
 - Quantitation of mammalian sialidase

REFERENCES

- Dharan, NJ et al (2009). *J Am Med Assoc* 301(10):1034-1041.
- Wetherill, NT et al (2003). *J Clin Microbiol* 41(2):742.
- Sheu, TG et al (2006). *Antimicrobial Agents & Chem* 52(9):3284.
- Cheng, PKC et al (2009). *Emerg Infect Dis* 15(6):966.
- Kawakami, C et al (2009). *Jpn J Infect Dis* 62:83.
- Matsuzaki, Y et al (2010). *Vir* 7:753.
- Leung, TWC et al (2005). *J Clin Virol* 46:298.
- An, J et al (2009). *J Med Chem* 52:2667.
- Hashem, AM et al (2009). *PLoS ONE* 4(12):e8350.
- Bowyer, NM et al (2008). *J Virol* 82(20):10052.
- Bright, R et al (2007). *Vaccine* 25:3871.
- NISN IC₅₀ Methods (<http://www.nisn.org/>):
 - Marathe, BM, and EA Govorkova (2010). St. Jude Children's Research Hospital.
 - Hurt, A. (2009). Standard Operating Procedure. WHO-Q25. WHO Collaborating Centre for Reference and Research on Influenza, Australia.
 - McKimm-Breschkin, JA. Enzyme inhibition assay for sensitivity to NA inhibitors. Fluorescence based assay (MUNANA).
 - Lackeyby, A. (2008). Determination of influenza virus susceptibility to neuraminidase inhibitors using a fluorescent substrate.
 - Jonges, M et al (2010). *J Clin Microbiol* 48(3):928-940.
 - Hassantouf, A et al (2010). *Vaccine* 28:790.

TRADEMARKS/LICENSING

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. Triton is a registered trademark of Dow Chemical Co. Prism is a registered trademark of GraphPad Software, Inc. For Research Use Only. Not intended for any animal or human therapeutic or diagnostic use.